

DESIGN

**V**eronique®  
c o m p r e s s i o n w e a r

Compression Garments Overview

Version 08.03.1.



### **Why should I use a medical compression garment?**

It is essential that the compression garment worn during the recovery period provides targeted compression, something that is not found in a standard girdle or athletic tights. Athletic tights, girdles and jog bras also do not have flat seams, an open crotch and are not intended for long-term wear.

Design Veronique® compression garments are very different from a pair of cycling shorts, a body girdle or bra you might see in your local store. Our garments are specifically designed for postoperative wear and have important features like flat, non-marking seams, inner panels that provide the targeted compression necessary during healing, and no labels or anything that will irritate your skin during extended wear.

### **About the company**

Design Veronique® was founded over 22 years ago by Veronica C. Smith and is committed to designing, manufacturing and providing the finest quality compression garments to doctors, hospitals and patients throughout the world. We will not compromise on product quality or service. All garments are designed from the finest quality components and materials and manufactured in our California headquarters to maintain the highest manufacturing and quality control standards.

### **Product line**

Design Veronique® has over 85 different styles of compression garments, available in three colors and a size range of XS through 4X for most styles. We also make custom garments if you don't find the style or size needed.

### **Finding the ideal garment for you**

We make a lot of garments, we have the ideal garment to help you recover from any procedure. Which garment is best for you? Consult our garment finder at: [www.designveronique.com/garmentfinder](http://www.designveronique.com/garmentfinder)

The garment finder works by having you answer three questions: The type of procedure you had, the exact procedure performed and the number of days post-operative you are. Based on your responses, the garment finder will suggest the correct garment style and either a zipper or non-zippered style.



### **Ordering & Questions**

Found the garment you need? You can order any of our products on our secure website: [www.DesignVeronique.com](http://www.DesignVeronique.com) or call 800.442.5800 and speak with one of our patient specialists. Please give us a call any time you have a question about one of our garments or a measurement.

# Features of Zippered, First Stage Garments

## *Designed for immediate postoperative use*

A zippered garment is worn during the first stage of recovery when you are swollen and tender. This phase typically ranges from immediately post-op to 10-14 days. The zippers allow you to change in and out of the garment with ease.

## *Sized to accommodate postoperative swelling and edema*

The zippered style sizes are based on your measurements before surgery. Although you will be smaller after your procedure, postoperative swelling typically means you'll remain at your pre-op size for the first few weeks after your surgery. Using pre-op measurements ensures a proper, comfortable fit.

We recommend that you have more than one garment. Immediately post-op the doctor will put you in a zippered garment, this garment might become stained or soiled. In order to maintain continuous compression, the second garment can be worn while the first garment is being laundered.



**#753**  
Zippered Above-Knee Girdle  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Lower abdominal, hip, and thigh procedures



**#85**  
Zippered Torso Brief  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal procedures



**#850**  
Zippered Abdominal Girdle  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal procedures



**#B85**  
Zippered High-Back Torso Brief with Bra  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal, back & breast procedures



**#B853**  
Zippered Above-Knee High-Back Girdle with Bra  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal, hip, thigh, back & breast procedures



**#855**  
Zippered Full-Body Girdle  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal, hip, thigh, and knee procedures

# Features of Non-Zippered, Second Stage Garments

## *Designed for body contouring and skin retracting*

As your swelling begins to subside, a non-zippered garment will provide more targeted support and aid in skin retraction during the final body contouring phase.

## *Targeted Compression Panels*

Non-zippered garments have strategically placed compression panels in key surgical areas to provide additional compression. The targeted compression provides more support and helps flush out any remaining fluids.

## *Designed for long-term wear*


Non-zippered garments are particularly comfortable for long-term wear and allow you to resume regular daily activities.



**#1453**  
Non-Zippered Above-Knee Girdle  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Lower abdominal, hip,  
and thigh procedures



**#165**  
Non-Zippered Torso Brief  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal procedures




**#1650**  
Non-Zippered Abdominal Girdle  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal procedures



**#1655**  
Non-Zippered Full-Body Girdle  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal, hip, thigh,  
and knee procedures



**#B165**  
Non-Zippered High-Back Torso  
Brief with Bra  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal, back & breast  
procedures



**#B1653**  
Non-Zippered Above-Knee High-  
Back Girdle with Bra  
Sizes: XS - 4X  
Colors: Black, White, Champagne  
Ideal for: Abdominal, hip, thigh,  
back & breast procedures

# Breast Wear

Visit [www.DesignVeronique.com](http://www.DesignVeronique.com) to see all Breast Wear garments.




**#450**  
D' Andrea — Optimum Support Cotton Knit Bra  
Sizes: XS – 3X  
Colors: Black, White  
Ideal for: Breast reduction and reconstruction



**#4510**  
Georgette—Fully Adjustable Implant Stabilizing Bra  
Sizes: XS - 3X  
Colors: Black, White, Champagne  
Ideal for: Breast augmentation, stabilizing  
and positioning implants



**#454**  
Allyssandra — 2" Band Cotton Knit Bra  
Sizes: XS – 3X  
Colors: Black, White, Champagne  
Ideal for: Breast augmentation and mastopexy



**#510-X**  
*Versatile Stabilizing Breast Wrap*  
Sizes: One Size  
Colors: White  
Ideal for: Breast augmentation, stabilizing,  
and positioning implants



**#457-Z**  
Yesmina-Z — Front Zippered Cotton Medical/Sports Bra  
Sizes: XS – 3X  
Colors: Black, White, Champagne  
Ideal for: Breast augmentation, stabilizing and position-  
ing implants



**#510-J**  
*3" Breast Wrap with Lace*  
Sizes: One Size  
Colors: Black, White  
Ideal for: Breast implant stabilization  
and positioning

## Upper Body Garments

Visit [www.DesignVeronique.com](http://www.DesignVeronique.com) to see all upper body garments.



**#B943**

Compression Arm Sleeve  
with Adjustable Cotton Knit Bra

Sizes: XS - 4X

Colors: Black, White, Champagne

Ideal for: Arm and breast procedures



**#BV933**

Long Arm Sleeve Bolero

Sizes: XS - 4X

Colors: Black, White

Ideal for: Arm, back, and breast procedures



**#350-9**

9" Double Panel Binder

Sizes: XS - 4X

Colors: White

Ideal for: Abdominal procedures

## Male Garments

Visit [www.DesignVeronique.com](http://www.DesignVeronique.com) to see all male garments.



**#640**

Adjustable Compression  
Vest

Sizes: S - 4X

Colors: Black, White

Ideal for: Chest procedures



**#1240**

Compression Vest

Sizes: S - 4X

Colors: Black, White

Ideal for: Chest procedures

## Facial Bands

Visit [www.DesignVeronique.com](http://www.DesignVeronique.com) to see all facial bands.



**#210**

Universal Facial Band<sup>®</sup>

Sizes: One Size

Colors: White, Champagne

Ideal for: Facial procedures



**#210-N**

Universal Facial and Neck Wrap

Sizes: M and L

Colors: White, Champagne

Ideal for: Facial and neck  
procedures